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English

dentona flexisplint®

General information regarding the work instructions

Caution storage:

- Liquid must be stored in a refrigerator at a temperature of 5 - 9°C.
- Powder must be stored at room temperature of 20 - 23°C.
- Both components should be cooled for approximately 15 minutes before mixing.

Please note:

- To prevent uncontrolled chemical reactions only use a mixing bowl and mixing spatula that do not come into contact with any other plastic. Repairs are performed in the customary manner.

The area is first roughened, then wetted with a little „dentona flexisplint®“ fluid and finally processed with fresh „dentona flexisplint®“ material.

Important:

- The mixing ratio and temperature must be strictly adhered to, as otherwise the plastic cannot cure correctly.

Care:

- Clean the dentona flexisplint® splint with a toothbrush and toothpaste under running water. Please do not allow to dry by evaporation! Disinfect with an alcohol-free cleaning agent before first use.

Pack pressing technique

- Prepare the model in the normal way, measure and mark the tooth equator. As dentona flexisplint® is a flexible material, the splint can be arranged slightly above the tooth equator.
- The splint to be created is waxed-up in the articulator taking into account the chewing movement (for example with pink wax sheet), and then bed the model into a cuvette. Now cool the required amount of powder in a refrigerator, (both components must be cooled before mixing).
- Isolate all exposed plaster parts of the embedding with alginate isolation. As soon as the isolation is dry the counterpart can be produced from the modelled wax splint. Once the counterpart wax has set completely hard the complete cuvette is warmed in a hot water bath. This aids opening the two halves of the cuvette.
- Remove the wax from the modelled wax splint, scald the model and the counterpart, water and isolate with alginate isolation. **TIP:** Severe undercut areas and interdental spaces should be sufficiently blocked off so that when the splint is completed it can be lifted without damage to / from the model.
- Remove the liquid from the refrigerator immediately prior to processing.
- Mixing ratio: 1 part powder / 0.75 parts fluid
IMPORTANT: Strictly adhere to this ratio!
 - 12 ml powder and 9 ml liquid are required for a small splint.
 - 16 ml powder and 12 ml liquid are required for a large splint.
- Measure the dentona flexisplint® material in the dosing cups (remove the liquid using the pipette provided), bring together in a separate rubber mixing cup, mix rapidly with the spatula excluding any bubbles for 5 to 10 seconds, cover and leave to stand.
- As soon as dentona flexisplint® adopts a paste-like consistency, (approximately 5 minutes), the plastic in the cuvette is placed over the row of teeth and the upper part of the cuvette is positioned. Any visible white spots in the stirred paste result from the production process and are no longer visible after polymerisation.
- Now quickly screw down under the press and then leave to polymerise in a bath of water at 2 bar and 50°C (maximum 53°C!).
- The cuvette must be completely cooled down before opening!

This is ensured by placing the cuvette in a cold bath of water for approximately 15 minutes.

- Now carefully lift the splint from the model and perform the usual finishing and polishing.

TIP: Before lifting off of the model and before setting in the patient’s mouth the splint must be made slightly flexible by exposing to lukewarm water.

Casting technique

- Prepare the model in the normal way, measure and mark the tooth equator. As dentona flexisplint® is a flexible material, the splint can be arranged slightly above the tooth equator.
- The splint to be created is waxed-up in the articulator taking into account the chewing movement (for example with pink wax sheet), and then bed the model into a cuvette. Now cool the required amount of powder in a refrigerator, (both components must be cooled before mixing).
- Now prepare a silicone counterpart (for example dentona® 1:1 gum) and allow to cure. Once the silicone counterpart is fully cured it can be carefully lifted off. Cut an excess channel free with a sharp scalpel on the tuber or both sides on the retromolar triangle.
- Remove the wax from the modelled wax splint, scald the model and the counterpart, water and isolate the still warm plaster model with alginate isolation (the silicone counterpart is then isolated with a thin layer of Vaseline). Once the isolation is dry and the model is completely cooled down, (simply place it in the refrigerator), place the counterpart back on the model and secure. **TIP:** Severe undercut areas and interdental spaces should be sufficiently blocked off so that when the splint is completed it can be lifted without damage to / from the model.
- Remove the liquid from the refrigerator immediately prior to processing.
- Mixing ratio: 1 part powder / 0.75 parts fluid
IMPORTANT: Strictly adhere to this ratio!
 - 12 ml powder and 9 ml liquid are required for a small splint.
 - 16 ml powder and 12 ml liquid are required for a large splint.
- Measure the dentona flexisplint® material in the dosing cups (remove the liquid using the pipette provided), bring together in a separate rubber mixing cup, mix rapidly with the spatula excluding any bubbles for 5 to 10 seconds, cover and leave to stand and immediately pour the free-flowing material into one side of the counterpart until the plastic can be seen on the other side of the counterpart.
- Now polymerise the model with the filled silicone counterpart as quickly as possible - for 45 minutes - in a pressure curing unit at 2 bar and a minimum of 50°C (maximum 53°C!).
- Then the workpiece is removed from the pressure curing unit and cooled in cold water before the finished splint is taken off.
- Now carefully lift the splint from the model and perform the usual finishing and polishing. **TIP:** Before lifting off of the model and before setting in the patient’s mouth the splint must be made slightly flexible by exposing to lukewarm water.

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