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Processing Instructions estheticast[®] exact

deutsch

english

français

italiano

español

estheticast® exact – Operating instructions

The estheticast® exact crown and bridge investment material is a carbon-free, phosphate-bonded precision investment material that can be used for both speed casting techniques as well as „conventional“ casting techniques. It is specially designed for precision castings in crown and bridge work with all common and usual alloys - it has been developed as high gold content, reduced gold content or base metal - as well as pressed ceramics.

estheticast® exact brief report:

- Suitable for „speed casting“ and controlled preheated, with mould and without mould ring.
- Extremely smooth castor pressed surfaces.
- Precisely fitting casting or pressing results even in large-span work such as „from ear to ear „or superstructures of implant restorations.
- Very easy investment due to its fine-grained consistency.
- Comfortable divestment.

Physical values

(Standard values at room temperature of 23°C and 100% liquid concentration)

ISO 15912:2006	Speed heating (Fitting after 25 minutes)	Controlled preheated (Fitting after 120 minutes)
Setting expansion	1,75 %	1,95 %
Thermal expansion	0,65 %	0,65 %
Total expansion (linear)	2,40 %	2,60 %
Setting time	~ 15 min.	
Compressive strength	8,5 MPa	
Processing time	6 min.	

Instructions for use:

1. Storage

Powder and liquid are stored at normal room temperature (23°C).

The material must be brought up to room temperature for at least 24 hours before use if it has been stored in temperatures that are too cold (such as a refrigerator). The storage and processing temperature of both components affect the setting time and the expansion of the investment.

INFO: Investment material liquid must be protected from temperatures of <5°C. Liquid that has been frozen cannot be used. Please be sure to obtain sufficient material before the onset of winter.

2. Processing temperature

Process the investment material at +/- 23°C (but at a minimum of 19°C).

INFO: Higher processing temperatures marginally reduce the processing time but improve the cast surface. Temperatures that are too low extend the processing time but can lead to roughness of the cast surface. Temperatures below 19°C extend the setting time.

3. Preparations prior to investment

For the waxing please use an insulating agent that does not leave an oily film.
Add the sprues following customary rules.

INFO: If you slightly relax the wax surface please ensure that the relaxant (e.g. coppie-flux) is completely dry before investment.

4. Mould ring sizes

Investment without a ring is also possible.

Line the metal mould with dento-fleece mould ring liner. Seal the fleece with a thin layer of Vaseline.

Open or wet fleece can change the liquid content of the mixture and therefore affect expansion.

1 mould = 1 liner dento fleece

3 mould = 1 liner dento fleece

6 mould = 2 liners dento fleece

9 mould = 2 liners dento fleece

The fleece must be flush with the top of the mould ring.

INFO: Take care with regard to mould size 1: Because of the small amount of material - due to the lower exothermic reaction - the cast may be slightly smaller. The best casting results are obtained with the mould size 3.

5. Ratio of powder / liquid to the mould size

Mould size	Powder (g)	Liquid (ml)
X3	150	39
X6	300	78
X9	450	117

INFO: Standard liquid ratio = 100 g / 26 ml

6. Expansion

The concentration of the liquid always depends on various factors

- Storage temperature of liquid and powder
- Mixing time
- Investment with or without pressure
- Amount of mixture

Mixing recommendation for optimum casting results*

Type of alloy	Muffle size	High gold content alloy > 70% Au		Reduced alloy > 50% Au		Pd – Base alloy		Base metal alloys		Pressed ceramics Total liquid 36 ml		General info	
		ml. Liquid	ml. dest. water	ml. Liquid	ml. dest. water	ml. Liquid	ml. dest. water	ml. Liquid	ml. dest. water	ml. Liquid	ml. dest. water	ml. Liquid	ml. dest. water
In/Onlays Pin structure Wax	3	19	20	20	19	26	13	30	9	19	17		
	6	38	40	36	42	52	26	60	18				
	9	57	60	54	63	78	39	90	37				
Crowns and bridges from wax	3	21	18	22	17	27	12	32	7	21	15		
	6	42	36	44	34	54	24	64	14				
	9	63	54	66	51	81	36	96	21				
Crowns and bridges, cone from pattern	3	23	16	25	14	29	10	35	4				
	6	46	32	50	28	58	20	70	8				
	9	69	48	75	42	87	30	105	12				
TK-outer parts from pattern	3	26	13	28	11	32	9	39 - 36					
	6	52	26	55	22	64	14	78 - 72					
	9	78	39	83	33	96	21	117 - 108					

If the expansion of the investment material for base materials is not sufficient, reduce the total amount of liquid to 36 ml / 150 g

*All information in this manual is based on intensive tests and a variety of casting trials. Various final results can be expected as each technician has a different way of working and different equipment. All values in this table are based on results of speed casting techniques. In conventional casting techniques, the mixing ratio may need to be adjusted slightly.

INFO: In order to increase the compressive strength, the press ceramic-mould is mixed with only 36 ml / 150 g. It is also recommended to prepare the press plunger from the „estheticast® exact“, so that the same expansion and compressive strength is achieved. To do so duplicate the original plunger and always fill the resulting shape simultaneously with the mould.

7. Processing and investment

Manually mix the appropriate amount of powder and liquid (see table) with a spatula until all the dry powder is wetted with liquid. Then place the mixing bowl in the vacuum mixer and stir for 60 s under vacuum. Proven mixing speeds are between 300 and 400 rev/min. Fill the mould with the investment material at a very low frequency. A probe or similar aid should be used to fill the crowns.

INFO: If you want to invest modeling resin (outer parts, bars, attachments etc.), you should post polymerise the plastic parts at 45 - 55°C prior to investment.

8. Setting down / Preheat

The mould must always be touched on the top. After which the mould can be placed on an undulating surface with the casting funnel pointing downwards.

- **Speed-cast:** After 25-30 min. – calculated from the first contact between powder/liquid – place in the preheated oven at 850°C for at least 45 minutes. Important: The muffle must still be hot! Now set the temperature to the desired final temperature. 850°C is sufficient for most alloys

NOTE: For ovens with bottom heating ensure sufficient distance (about 1 cm) between the mould and the bottom

After 45 min. in a hot oven (No. 9 mould 60 min., no. 1 mould 35 min.) swiftly pour off.

- **Conventional preheating:** After the investment material is set, place the mould into a cold oven and heat at a rate of temperature increase of 3-5°C / min. Allow for a holding time of 30 min at 290°C and 590°C. Hold the final temperature (for alloy) for at least 45 min. (No. 9 mould 60 min., no. 1 mould 35 min.).
- **Special feature of modelling resin:** After 25 mins place the invested mould into an oven preheated to 360°C for 30 mins in order to soften the modelling resin. Only then increase the heat up to the final temperature.

9. Casting

After the recommended preheating time of 45 min. cast in the usual way. The preheat period is increased by 10 min per mould if more than one mould is preheated in the oven.

10. Cool

Let the mould cool slowly at room temperature. NEVER chill.

The above information is correct to the best of our knowledge and careful consideration. It reflects the present state of the art. We guarantee the faultless quality of our products, but are not liable for further processing results, which usually arise outside our sphere of influence.

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dentona AG | Otto-Hahn-Straße 27 | 44227 Dortmund | Germany