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ACRYLIC RESIN FOR TEMPORARY CROWNS AND BRIDGES

FOR DENTAL USE ONLY

Composition

Self-curing acrylic resin consisting of:

Powder: polymethylmethacrylate

Liquid: methyl methacrylate, dimethacrylate

Intended use

Acrylic resin for temporary crowns and bridges, immediate crowns brush-modeling, for the gypsum or silicone key pouring technique, for rapid repairing of acrylic crowns and bridges.

GMDN Code 31783

RESPAL NF complies with reference technical standards ISO 10477:2004 - Tipo 1.

Warnings

- ◆ Methacrylate-based resins may cause sensitization in contact with skin and mucous membranes. It is not recommended for use in the event of a suspected allergy to the components. Check the patient's reaction during the first period of exposure to the product. Not polymerize in the mouth or in contact with oral tissues.
- ◆ The methyl methacrylate contained in the monomer is highly flammable (flash point +9°C). Store and use away from sources of heat and ignition. Do not smoke during use.
- ◆ Avoid inhaling the vapors of the monomer because they can cause irritation to eyes and respiratory tract.
- ◆ Avoid prolonged contact with uncured product that, in susceptible subjects, may cause sensitization and irritation of the skin.
- ◆ Read and carefully follow the safety instructions contained in the *Product Safety Data Sheets*, which are available on our web site at www.spd-dental.com

Available colours

TAKILON BB is available in 7 different colours selected to give the product a natural aesthetic appearance. Enamels, that are available separately, provide an added possibility of customizing the colors as needed.

- | | | | |
|------------------------------------|----------------------------------|--------------------------------|-----------------------------------|
| ◆ Light values:
Colours 1, 2, 3 | ◆ Medium values:
Colours 4, 5 | ◆ Dark values:
Colours 6, 7 | ◆ Enamels: Blue
(A), White (B) |
|------------------------------------|----------------------------------|--------------------------------|-----------------------------------|

Preparing the model

Isolate the model with a gypsum-resin alginic insulator.

Dosing

Ideal mixing proportion: 1 ml of liquid every 2,5 g of powder.

Empirical method: 3 parts powder to 1 part liquid by volume.

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[8]

Preparing the mixture

Pour the liquid and the powder components into a porcelain or glass mortar in the correct dosage. Mix thoroughly with a stainless steel spatula for about 30 seconds.

TAKILON BB is ready to use as soon as it separates easily from the spatula.

Application

Apply the mixture directly on the model with the usual laboratory techniques for the construction of temporary crowns and bridges.

Working time

TAKILON BB is self-curing, the complete hardening of the mass at room temperature of 23°C will occur within a median time of 8 minutes.

Curing

To obtain a highly homogeneous final product, with great firmness and strength, immerse the model in **POL-01 Pressure pot** for 20 minutes at a temperature of 40°C and a pressure of about 2 bar. The model should be completely immersed in water. The polymerization in the open air, although achievable, is strongly not recommended.

Finishing

Once the polymerization is carried out, remove the artifact from the model. Finish and polish the artifact in the ordinary way, by using burs, separators, wheels and polishing pastes recommended by the manufacturer for use with acrylics. Do not over-finishing and polishing: the heat generated during such operations can damage the material. The use of the same tools with different materials and colors can alter the color of the product.

Storage

Store the product in cool and ventilated room, protected from light and away from heat sources. Keep out of reach of children and unauthorized personnel.

Do not litter and do not spill the monomer into drains. Dispose of in accordance with local and national regulations.

Sales formats

Kit containing 7 Powders 100 g + 2 Liquids 250 ml

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